



Patrol Pump Technical Sheet

A great solution for the
non-ambulatory tube-fed patient.



EASY-TO-USE

SAFE

Simple Set up
Free-flow Alarms
Programmable HOLD Timer
Lock-out Feature

- ✓ Free flow prevented with the **Safe-T-Valve™**
- ✓ **System self check**
- ✓ Useful **1-to-90 minute programmable HOLD**
- ✓ Tampering prevented with the **lockout feature**
- ✓ Reduced risk of error with the **memory feature**
- ✓ **8-hour battery operation**
- ✓ Important **safety alarms** such as select rate, free flow, dose fed, low battery and failed self test
- ✓ Convenient **special holder** for the feeding set covers
- ✓ User-friendly **control dial**
- ✓ Easy-to-read **large LED display**
- ✓ Easy to **set up** and **prime**
- ✓ **Easy to use** and **easy to teach** staff and patients
- ✓ Suitable for adult or pediatric patients

Patrol Pump Sets

CLOSED SYSTEM

Formula is prefilled into feeding container. Reduced handling decreases risk of microbial contamination.
 Ross Ready-to-Hang® containers can hang up to 48 hours when one feeding set and clean techniques are used.

PUMP SETS

Product Code	Units/Case	
52040-U30	30	Piercing Pin



Ross Ready-to-Hang® containers



OPEN SYSTEM

Formula needs to be decanted from the original container into the feeding container.

- Ready-to-use formula can hang up to 8-12 hours.
- Reconstituted formula or formula with additives such as modular components can hang up to 4 hours.
- All pump sets for an open system are to be changed within 24 hours.

PUMP SETS

OPTION 1

NEW



Product Code	Units/Case	
52048-030	30	Easy-Feed® preattached pump set

OPTION 2



Product Code	Units/Case	
52042-U30	30	Top-Fill® with preattached pump set

OPTION 3



Product Code	Units/Case	
52040-U30	30	Piercing Pin



Product Code	Units/Case	
90-024	24	Top-Fill®